

[Header]							
IEMFileVersion	4						
Investigator Name	Sample						
Experiment Name	Sample						
Date	00/00/0000						
Workflow	GenerateFASTQ						
Application	FASTQ Only						
Assay	TruSeq LT						
Description	Sample						
Chemistry	Default						
[Reads]							
	151						
	151						
[Settings]							
ReverseComplement	0						
Adapter	AGATCGGAAGAGCACACGTCTGAACTCCAGTCA						
AdapterRead2	AGATCGGAAGAGCGTCGTGTAGGGAAAGAGTGT						
[Data]							
Sample_ID	Sample_Name	Sample_Plate	Sample_Well	I7_Index_ID	index	Sample_Proje	Description
Sample Name				Adapter 1	AACGTG		
Sample Name				Adapter 2	CAGATC		
Sample Name				Adapter 3	GTACGC		
Sample Name				Adapter 4	TATCAG		
Sample Name				Adapter 5	TCTTCA		
Sample Name				Adapter 6	CTAAGG		
Sample Name				Adapter 7	GAGTTA		
Sample Name				Adapter 8	AAGGTA		
Sample Name				Adapter 9	ACGCTC		
Sample Name				Adapter 10	AGATCG		
Sample Name				Adapter 11	ATCCTG		
Sample Name				Adapter 12	GACTAG		
Sample Name				Adapter 13	ATGCCT		
Sample Name				Adapter 14	CATCAA		
Sample Name				Adapter 15	AGTACA		
Sample Name				Adapter 16	GAATCT		
Sample Name				Adapter 17	GATAGA		
Sample Name				Adapter 18	GCTCGG		
Sample Name				Adapter 19	GTCTGT		
Sample Name				Adapter 20	TCCGTC		
Sample Name				Adapter 21	TGAAGA		
Sample Name				Adapter 22	AACGCT		
Sample Name				Adapter 23	ACGTAT		
Sample Name				Adapter 24	AGTCAC		