

<b>Adapter ID</b>	<b>6bp Index Sequence</b>
Adapter 1	<b>AACGTG</b>
Adapter 2	<b>CAGATC</b>
Adapter 3	<b>GTACGC</b>
Adapter 4	<b>TATCAG</b>
Adapter 5	<b>TCTTCA</b>
Adapter 6	<b>CTAAGG</b>
Adapter 7	<b>GAGTTA</b>
Adapter 8	<b>AAGGTA</b>
Adapter 9	<b>ACGCTC</b>
Adapter 10	<b>AGATCG</b>
Adapter 11	<b>ATCCTG</b>
Adapter 12	<b>GACTAG</b>
Adapter 13	<b>ATGCCT</b>
Adapter 14	<b>CATCAA</b>
Adapter 15	<b>AGTACA</b>
Adapter 16	<b>GAATCT</b>
Adapter 17	<b>GATAGA</b>
Adapter 18	<b>GCTCGG</b>
Adapter 19	<b>GTCTGT</b>
Adapter 20	<b>TCCGTC</b>
Adapter 21	<b>TGAAGA</b>
Adapter 22	<b>AACGCT</b>
Adapter 23	<b>ACGTAT</b>
Adapter 24	<b>AGTCAC</b>