

Seraseq[®] ctDNA Complete[™]

Seraseq ctDNA Complete reference materials are purpose-built to accelerate the thorough validation and clinical implementation of NGS-based ctDNA assays. Designed with the most complete coverage of actionable variant types in a single reference sample, they're widely compatible with a variety of commercially available NGS-based liquid biopsy assays.

NGS-Based Assay Compatibility					
Archer [®] Reveal ctDNA [™] 28 Kit	Oncomine [®] Pan-Cancer Cell-Free Assay	Oncomine Lung cfDNA Assay	Oncomine Colon cfDNA	UltraSEEK [™] Colon Panel	UltraSEEK Lung Panel
AKT1 ●	AKT1 ●	ALK ●▲	AKT1 ●	BRAF ●	BRAF ●
ALK ●▲	ALK ●▲	BRAF ●	BRAF ●	EGFR ●◇	EGFR ●◇
BRAF ●	BRAF ●	EGFR ●◇	EGFR ●◇	KRAS ●	ERBB2 ◇■
EGFR ●◇	EGFR ●◇	ERBB2 ◇■	ERBB2 ◇■	NRAS ●	KRAS ●
ERBB2 ◇■	ERBB2 ◇■	KRAS ●	KRAS ●	PIK3CA ●◇	PIK3CA ●◇
KIT ●	KIT ●	MET ■	NRAS ●		
KRAS ●	KRAS ●	NRAS ●	PIK3CA ●◇		
MET ■	MET ■	PIK3CA ●◇			
NRAS ●	MYC ■	ROS1 ▲			
PIK3CA ●■	NRAS ●				
RET ▲	PIK3CA ●◇				

Key	
SNV	●
INDEL	◇
CNV	■
Fusion	▲

Learn more at www.seracare.com/ctDNA